





CDKAL1 Antibody

Product Code	CSB-PA716574LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VV42
Immunogen	Recombinant Human Threonylcarbamoyladenosine tRNA methylthiotransferase protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Threonylcarbamoyladenosine tRNA methylthiotransferase (EC 2.8.4.5) (CDK5 regulatory subunit-associated protein 1-like 1) (tRNA-t(6)A37 methylthiotransferase), CDKAL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDKAL1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA716574LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA716574LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)