



Fbxw10 Antibody, Biotin conjugated

Product Code	CSB-PA725549LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SUS0
Immunogen	Recombinant Mouse F-box/WD repeat-containing protein 10 protein (701-1030AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 10 (F-box and WD-40 domain-containing protein 10), Fbxw10
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Fbxw10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.