



Gfer Antibody

Product Code	CSB-PA733907HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63042
Immunogen	Recombinant Rat FAD-linked sulfhydryl oxidase ALR protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAD-linked sulfhydryl oxidase ALR (EC 1.8.3.2) (Augmenter of liver regeneration), Gfer, Alr
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Gfer
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.