



SLC16A9 Antibody, HRP conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA742402LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7RTY1 |
| Immunogen | Recombinant Human Monocarboxylate transporter 9 protein (185-304AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Monocarboxylate transporter 9 (MCT 9) (Solute carrier family 16 member 9), SLC16A9, C10orf36 MCT9 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SLC16A9 |