



SLC16A9 Antibody, FITC conjugated

Product Code	CSB-PA742402LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTY1
Immunogen	Recombinant Human Monocarboxylate transporter 9 protein (185-304AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 9 (MCT 9) (Solute carrier family 16 member 9), SLC16A9, C10orf36 MCT9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A9