



SLC16A9 Antibody, Biotin conjugated

Product Code	CSB-PA742402LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTY1
Immunogen	Recombinant Human Monocarboxylate transporter 9 protein (185-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 9 (MCT 9) (Solute carrier family 16 member 9), SLC16A9, C10orf36 MCT9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.