



# WAPL Antibody, Biotin conjugated

Product Code	CSB-PA772025LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5K2
Immunogen	Recombinant Human Wings apart-like protein homolog protein (1037-1190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wings apart-like protein homolog (Friend of EBNA2 protein) (WAPL cohesin release factor), WAPL, FOE KIAA0261 WAPAL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WAPL