



# SLC15A3 Antibody, Biotin conjugated

Product Code	CSB-PA815578LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY34
Immunogen	Recombinant Human Solute carrier family 15 member 3 protein (253-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 15 member 3 (Osteoclast transporter) (Peptide transporter 3) (Peptide/histidine transporter 2), SLC15A3, OCTP PHT2 PTR3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC15A3