



WDFY1 Antibody, HRP conjugated

Product Code	CSB-PA816890LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWB7
Immunogen	Recombinant Human WD repeat and FYVE domain-containing protein 1 protein (261-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat and FYVE domain-containing protein 1 (FENS-1) (Phosphoinositide-binding protein 1) (WD40- and FYVE domain-containing protein 1) (Zinc finger FYVE domain-containing protein 17), WDFY1, KIAA1435 WDF1 ZFYVE17
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	WDFY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.