



TP53INP1 Antibody, HRP conjugated

Product Code	CSB-PA822174LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A56
Immunogen	Recombinant Human Tumor protein p53-inducible nuclear protein 1 protein (138-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible nuclear protein 1 (Stress-induced protein) (p53-dependent damage-inducible nuclear protein 1) (p53DINP1), TP53INP1, P53DINP1 SIP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53INP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.