



TP53INP1 Antibody, HRP conjugated

Product Code	CSB-PA822174LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A56
Immunogen	Recombinant Human Tumor protein p53-inducible nuclear protein 1 protein (138-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible nuclear protein 1 (Stress-induced protein) (p53-dependent damage-inducible nuclear protein 1) (p53DINP1), TP53INP1, P53DINP1 SIP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53INP1