



BRSK1 Antibody, HRP conjugated

Product Code	CSB-PA837438LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDC3
Immunogen	Recombinant Human Serine/threonine-protein kinase BRSK1 protein (426-591AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase BRSK1 (EC 2.7.11.1) (Brain-selective kinase 1) (EC 2.7.11.26) (Brain-specific serine/threonine-protein kinase 1) (BR serine/threonine-protein kinase 1) (Serine/threonine-protein kinase SAD-B) (Synapses of Amphids Defective homolog 1) (SAD1 homolog) (hSAD1), BRSK1, KIAA1811 SAD1 SADB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BRSK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.