



PNKD Antibody, Biotin conjugated

Product Code	CSB-PA843154LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N490
Immunogen	Recombinant Human Probable hydrolase PNKD protein (78-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable hydrolase PNKD (EC 3.-.-.) (Myofibrillogenesis regulator 1) (MR-1) (Paroxysmal nonkinesinogenic dyskinesia protein) (Trans-activated by hepatitis C virus core protein 2), PNKD, KIAA1184 MR1 TAHCCP2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PNKD