



PNKD Antibody, Biotin conjugated

Product Code	CSB-PA843154LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N490
Immunogen	Recombinant Human Probable hydrolase PNKD protein (78-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable hydrolase PNKD (EC 3.-.-.) (Myofibrillogenesis regulator 1) (MR-1) (Paroxysmal nonkinesinogenic dyskinesia protein) (Trans-activated by hepatitis C virus core protein 2), PNKD, KIAA1184 MR1 TAHCCP2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PNKD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.