



KCNG4 Antibody, Biotin conjugated

Product Code	CSB-PA844728LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDN1
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily G member 4 protein (1-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily G member 4 (Voltage-gated potassium channel subunit Kv6.4), KCNG4, KCNG3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.