



KCNG4 Antibody, Biotin conjugated

Product Code	CSB-PA844728LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDN1
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily G member 4 protein (1-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily G member 4 (Voltage-gated potassium channel subunit Kv6.4), KCNG4, KCNG3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNG4