



GRHL3 Antibody, Biotin conjugated

Product Code	CSB-PA844730LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE85
Immunogen	Recombinant Human Grainyhead-like protein 3 homolog protein (31-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Grainyhead-like protein 3 homolog (Sister of mammalian grainyhead) (Transcription factor CP2-like 4), GRHL3, SOM TFCP2L4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GRHL3