





PPP1R13B Antibody

Product Code	CSB-PA853449LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KQ4
Immunogen	Recombinant Human Apoptosis-stimulating of p53 protein 1 protein (375-675AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis-stimulating of p53 protein 1 (Protein phosphatase 1 regulatory subunit 13B), PPP1R13B, ASPP1 KIAA0771
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PPP1R13B

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA853449LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA853449LA01HU at dilution of 1:100

