



GNPDA2 Antibody, FITC conjugated

Product Code	CSB-PA855058LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDQ7
Immunogen	Recombinant Human Glucosamine-6-phosphate isomerase 2 protein (28-154AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucosamine-6-phosphate isomerase 2 (EC 3.5.99.6) (Glucosamine-6-phosphate deaminase 2) (GNPDA 2) (GlcN6P deaminase 2) (Glucosamine-6-phosphate isomerase SB52), GNPDA2, GNP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNPDA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.