



Grk2 Antibody

Product Code	CSB-PA858153LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99MK8
Immunogen	Recombinant Mouse Beta-adrenergic receptor kinase 1 protein (1-689AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein-coupled receptor kinase 2), Grk2, Adrbk1
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Grk2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.