




MYD88 Antibody

Product Code	CSB-PA859945LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99836
Immunogen	Recombinant Human Myeloid differentiation primary response protein MyD88 protein (1-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloid differentiation primary response protein MyD88, MYD88
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYD88

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA859945LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)