




DHX36 Antibody

Product Code	CSB-PA862002LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2U1
Immunogen	Recombinant Human ATP-dependent RNA helicase DHX36 protein (68-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase DHX36 (EC 3.6.4.12) (EC 3.6.4.13) (DEAH box protein 36) (G4-resolvase 1) (G4R1) (MLE-like protein 1) (RNA helicase associated with AU-rich element ARE), DHX36, DDX36 KIAA1488 MLEL1 RHAU
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DHX36

Image



Immunofluorescent analysis of U251 cells using CSB-PA862002LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.