



S1PR5 Antibody, HRP conjugated

Product Code	CSB-PA867132LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H228
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 5 protein (309-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingosine 1-phosphate receptor 5 (S1P receptor 5) (S1P5) (Endothelial differentiation G-protein-coupled receptor 8) (Sphingosine 1-phosphate receptor Edg-8) (S1P receptor Edg-8), S1PR5, EDG8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S1PR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.