



SLCO1C1 Antibody, FITC conjugated

Product Code	CSB-PA868355LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYB5
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1C1 protein (297-348AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1C1 (Organic anion transporter F) (OATP-F) (Organic anion transporter polypeptide-related protein 5) (OAT-RP-5) (OATPRP5) (Organic anion-transporting polypeptide 14) (OATP-14) (Solute carrier family 21 member 14) (Thyroxine transporter), SLCO1C1, OATP14 OATP1C1 OATPF SLC21A14
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLCO1C1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.