



MAN2C1 Antibody, HRP conjugated

Product Code	CSB-PA878885LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTJ4
Immunogen	Recombinant Human Alpha-mannosidase 2C1 protein (746-969AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-mannosidase 2C1 (EC 3.2.1.24) (Alpha mannosidase 6A8B) (Alpha-D-mannoside mannosidase) (Mannosidase alpha class 2C member 1), MAN2C1, MANA MANA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAN2C1