



MAN2C1 Antibody, HRP conjugated

Product Code	CSB-PA878885LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTJ4
Immunogen	Recombinant Human Alpha-mannosidase 2C1 protein (746-969AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-mannosidase 2C1 (EC 3.2.1.24) (Alpha mannosidase 6A8B) (Alpha-D-mannoside mannosidase) (Mannosidase alpha class 2C member 1), MAN2C1, MANA MANA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAN2C1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.