





FUZ Antibody

Product Code	CSB-PA883593LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT04
Immunogen	Recombinant Human Protein fuzzy homolog protein (280-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein fuzzy homolog, FUZ, FY
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FUZ

Image



Immunohistochemistry analysis of human melanoma using CSB-PA883593LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA883593LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.