




DIAPH3 Antibody

Product Code	CSB-PA885723LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSV4
Immunogen	Recombinant Human Protein diaphanous homolog 3 protein (642-849AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein diaphanous homolog 3 (Diaphanous-related formin-3) (DRF3) (MDia2), DIAPH3, DIAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DIAPH3

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA885723LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)