



DTL Antibody, HRP conjugated

Product Code	CSB-PA889160LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZJ0
Immunogen	Recombinant Human Denticleless protein homolog protein (393-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Denticleless protein homolog (DDB1- and CUL4-associated factor 2) (Lethal(2) denticleless protein homolog) (Retinoic acid-regulated nuclear matrix-associated protein), DTL, CDT2 CDW1 DCAF2 L2DTL RAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DTL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.