



DTL Antibody, Biotin conjugated

Product Code	CSB-PA889160LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZJ0
Immunogen	Recombinant Human Denticleless protein homolog protein (393-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Denticleless protein homolog (DDB1- and CUL4-associated factor 2) (Lethal(2) denticleless protein homolog) (Retinoic acid-regulated nuclear matrix-associated protein), DTL, CDT2 CDW1 DCAF2 L2DTL RAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DTL