



PCMP-H63 Antibody, Biotin conjugated

Product Code	CSB-PA890450YD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SUU7
Immunogen	Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At4g32450, mitochondrial protein (111-537AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pentatricopeptide repeat-containing protein At4g32450, mitochondrial, PCMP-H63
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	PCMP-H63
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.