



TBX21 Antibody, HRP conjugated

Product Code	CSB-PA890744LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL17
Immunogen	Recombinant Human T-box transcription factor TBX21 protein (82-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-box transcription factor TBX21 (T-box protein 21) (T-cell-specific T-box transcription factor T-bet) (Transcription factor TBLYM), TBX21, TBET TBLYM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TBX21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.