




SNX9 Antibody

Product Code	CSB-PA896917LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X1
Immunogen	Recombinant Human Sorting nexin-9 protein (85-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-9 (SH3 and PX domain-containing protein 1) (Protein SDP1) (SH3 and PX domain-containing protein 3A), SNX9, SH3PX1 SH3PXD3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNX9

Image



Western Blot

Positive WB detected in: Hela whole cell lysate, K562 whole cell lysate


All lanes: SNX9 antibody at 6.2µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67 kDa


Observed band size: 67 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA896917LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA896917LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in: HeLa whole cell lysate, K562 whole cell lysate, A549 whole cell lysate, U87 whole cell lysate, U251 whole cell lysate, 293 whole cell lysate

All lanes: SNX9 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67 kDa

Observed band size: 67 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.