



SNX9 Antibody, HRP conjugated

Product Code	CSB-PA896917LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X1
Immunogen	Recombinant Human Sorting nexin-9 protein (85-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-9 (SH3 and PX domain-containing protein 1) (Protein SDP1) (SH3 and PX domain-containing protein 3A), SNX9, SH3PX1 SH3PXD3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNX9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.