



SNX9 Antibody, Biotin conjugated

Product Code	CSB-PA896917LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X1
Immunogen	Recombinant Human Sorting nexin-9 protein (85-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-9 (SH3 and PX domain-containing protein 1) (Protein SDP1) (SH3 and PX domain-containing protein 3A), SNX9, SH3PX1 SH3PXD3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNX9