



HRH3 Antibody, HRP conjugated

Product Code	CSB-PA897568LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5N1
Immunogen	Recombinant Human Histamine H3 receptor protein (236-332AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histamine H3 receptor (H3R) (HH3R) (G-protein coupled receptor 97), HRH3, GPCR97
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HRH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.