



Recombinant Human Frizzled-7 (FZD7), partial

Product Code	CSB-BP009110HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75084
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QPYHGEKGISVPDHGFCQPISIP LCTDIAYNQ TILPNLLGHTNQEDAGLEVHQF YPLVKV QCSPELRFFLC SMYAPVCTVLDQAIPPCRSLCERARQGCEALMNKFGFQWPE RLRCENFP VHGAGEICV GQNTSDGSGGPGGGPTAYPTAPYLPDLPFTALPPGASDGRGRP AFPFSCPR QLKVPPYLG YRFLGERDCGAPCEPGRANGLMYFKEEERFARLW
Source	Baculovirus
Target Names	FZD7
Protein Names	Recommended name: Frizzled-7 Short name= Fz-7 Short name= hFz7 Alternative name(s): FzE3
Expression Region	33-256aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial, the first extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.