




CACNA2D3 Antibody

Product Code	CSB-PA741101
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZS8
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen Species	Homo sapiens (Human)
Target Names	CACNA2D3

Image



Western blot analysis of 1) Mouse Brain Tissue, 2) Human Brain Tissue, with CaV α 2 δ 3 Rabbit pAb diluted at 1:2,000.