





CACNA1C Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA775374 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13936 |
| Immunogen | Synthetic Peptide |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | IHC,ELISA;IHC:1:100-200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen. |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | CACNA1C |

Image



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Cav pan α 1Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Cav pan α 1Rabbit pAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.