




CCKBR Antibody

Product Code	CSB-PA924152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32239
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Alias	CCKBR; CCKRB; Gastrin/cholecystokinin type B receptor; CCK-B receptor; CCK-BR; Cholecystokinin-2 receptor; CCK2-R
Immunogen Species	Homo sapiens (Human)
Target Names	CCKBR

Image



Western blot analysis of 1) Mouse Brain Tissue, 2)Rat Brain Tissue, 3)Hela whole cell lysate, with CCKBR Rabbit pAb diluted at 1:2,000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.