




GALR2 Antibody

Product Code	CSB-PA627841
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43603
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Alias	GALR2; GALNR2; Galanin receptor type 2; GAL2-R; GALR-2
Immunogen Species	Homo sapiens (Human)
Target Names	GALR2

Image



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue with GALR2 Rabbit pAb diluted at 1:2,000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.