





KCNK9 Antibody

Product Code	CSB-PA836402
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPC2
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:1000-2000, IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen Species	Homo sapiens (Human)
Target Names	KCNK9

Image



Western blot analysis of 1) Mouse BrainTissue, 2)Rat Brain Tissue with KCNK9 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using KCNK9 (TASK-3) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using KCNK9 (TASK-3) Rabbit pAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.