




# KCNA5 Antibody

<b>Product Code</b>	CSB-PA914616
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P22460
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB,ELISA;WB:1:1000-2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Alias</b>	KCNA5; Potassium voltage-gated channel subfamily A member 5; HPCN1; Voltage-gated potassium channel HK2; Voltage-gated potassium channel subunit Kv1.5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	KCNA5

## Image



Western blot analysis of 293T with KV1.5 Rabbit pAb diluted at 1:2,000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.