





KCNA10 Antibody

Product Code	CSB-PA988090
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16322
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:1000-2000, IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen Species	Homo sapiens (Human)
Target Names	KCNA10


Image



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue with Kv1.8 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Kv1.8 Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Kv1.8 Rabbit pAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.