





KCNV2 Antibody

Product Code	CSB-PA681990
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDN2
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:1000-2000, IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Alias	KCNV2; Potassium voltage-gated channel subfamily V member 2; Voltage-gated potassium channel subunit Kv8.2
Immunogen Species	Homo sapiens (Human)
Target Names	KCNV2

Image



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue with Kv11.1 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Kv11.1 Rabbit pAb diluted at 1:200.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.