




TOP3A Antibody

Product Code	CSB-PA698044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13472
Immunogen	Synthesized peptide derived from the Internal region of Human Topo III α . at AA range: 350-430
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	TOP3A; TOP3; DNA topoisomerase 3-alpha; DNA topoisomerase III alpha
Immunogen Species	Homo sapiens (Human)
Target Names	TOP3A

Image



Western blot analysis of 293T using Topo III α antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.