




ANKRD20A1 Antibody

Product Code	CSB-PA865016
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TYW2
Immunogen	Synthesized peptide derived from the Internal region of Human A20.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD20A1

Image



Western blot analysis of 3T3 using A20 antibody.. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.