




AXL Antibody

Product Code	CSB-PA297707
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30530
Immunogen	Synthesized peptide derived from Human Axl around the non-phosphorylation site of Y691.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Immunogen Species	Homo sapiens (Human)
Target Names	AXL

Image



Western blot analysis of 293T using Axl antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.