




FGFR4 Antibody

Product Code	CSB-PA604644
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22455
Immunogen	Synthesized peptide derived from the Internal region of Human FGFR-4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR4

Image



Western blot analysis of 293T 3T3 lysis using FGFR4 antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000