




MME Antibody

Product Code	CSB-PA897259
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08473
Immunogen	Recombinant Protein of CD10 of MME
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC,ELISA;IHC:1:100-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	MME
Immunogen Species	Homo sapiens (Human)
Target Names	MME

Image



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using CD10 Rabbit pAb diluted at 1:500.